

DUCTLESS FUME HOOD. CLASSIC RANGE

# CRUMAECO<sup>2</sup>



**ECO<sup>2</sup>** is the fumehood that Cruma has developed responding to the demands of our clients.

Our R&D department has created this new model following two main premises: **Economy** and **Ecology**. Therefore its name. **ECO<sup>2</sup>**, meeting the needs of a market in evolution and constantly demanding new products, respects both the environment and the lab user with the same performance, functionality and quality of all Cruma products.

As a matter of fact, ECO<sup>2</sup> has passed the same quality and security tests of all other Cruma hoods, satisfying the requirements of all clients.




## ESPECIFICACIONES TÉCNICAS

Number of filtration columns	1
Number of filters	1 to 2
Number of IP44 fans	1
Average volume of treated air	160 m <sup>3</sup> /h
Average face velocity	0,50 m/s
Internal volume of the cabinet	0,287 m <sup>3</sup>
Renewals inside the cabinet / min	9,2
Total electrical power consumption	60 W
Voltage-Frequency	110-220 V - 50-60 Hz
Light intensity	1500 lm
Noise level	55 dB
<b>Packaging:</b> 100% recycled wooden box with international phytosanitary certificate	Volume 0,42 m <sup>3</sup> Weight 77 Kg

## SIZES (MM)

External			Internal		
Width	Depth	Height	Width	Depth	Height
798	600	1012	778	569	760

Do you have any *questions* or need *technical assistance*?



Contact your distributor or call us if you have any questions or need technical support, spare parts, maintenance service... ☎ +34 93 370 61 62

## SIZES (mm)

